Instructions for Use

 Manufactured by:

Z-Medica®, LLC

4 Fairfield Boulevard

Wallingford, CT 06492 USA

Tel: 1.877.750.0504, +1.203.294.0000

Fax: 1.800.343.8656, +1.203.303.7216

QuikClot.com

Patents: Z-Medica.com/patents.aspx ©

2021 Z-Medica, LLC

All rights reserved.

101683-03-WEB Rev A 5/11/2021

2797

EC REP

Emergo Europe
Prinsessegracht 20
2514 AP

The Hague
The Netherlands

Rx ONLY

REF 374

Please contact Z-Medica Customer Service to determine if this product is available for sale in your market / country.

Table of Contents

English.....	3-4
Preclinical Data.....	5
Clinical Data.....	6-8
Czech.....	9-10
French.....	11-12
German.....	13-14
Greek.....	15-16
Italian.....	17-18
Polish.....	19-20
Portuguese.....	21-22
Russian.....	23-24
Spanish.....	25-26
Swedish.....	27-28
Turkish.....	29-30

Description

QuikClot® Radial® consists of a roll dressing and an adhesive bandage. The roll is a soft, white, sterile, hydrophilic hemostatic dressing. This product contains kaolin.

Indications

QuikClot Radial is applied topically as an adjunct to manual compression and is indicated for the local management and control of surface bleeding from vascular access sites, percutaneous catheters or tubes utilizing introducer sheaths up to 12 Fr. or up to 7 Fr. for patients on drug/induced anti-coagulation treatment.

Contraindications

QuikClot Radial has not been tested on patients with bleeding disorders due to underlying disease (liver, kidney, or others) and is not indicated for these populations.

Warnings

QuikClot Radial has been tested in clinical trials and its efficacy has been shown only in patients treated with anticoagulation medications: heparin, clopidogrel bisulfate and warfarin. The efficacy of QuikClot Radial in the presence of other anticoagulation medications is not known.

Precautions

- For external use only
- Sterility not guaranteed if package is damaged or opened
- Discard if package is damaged
- Avoid contact with eyes
- Use aseptic techniques
- Do not re-sterilize
- Reuse will cause risk of infection and loss of efficacy
- Longer compression time may be required for patients who are hypertensive, obese, or on drug-induced anticoagulant therapies
- The adhesive bandage is not intended for use on fragile or compromised skin unless used in conjunction with protective films
- Do not use when an individual is sensitive to any of the components in the adhesive bandage
- In patients with either very large or very small wrists, where adequate pressure cannot be achieved by the adhesive bandage, other means such as medical tape or manual compression should be used to maintain firm pressure on the dressing
- If insufficient pressure is applied over the puncture site, bleeding may not be controlled

Storage Conditions


Keep dry. Keep from heat (including storage in direct sunlight or in direct contact with heat sources). Store at temperatures of 77°F (25°C) or less.

Rx ONLY means prescription use only


SBS means sterile barrier system

MD means medical device

Instructions for Use


(1) Remove hemostatic roll from package. Place roll into sterile field using aseptic technique. Do not wet the roll with saline before using it.



(2) Remove the adhesive bandage from package. Grasp the two edges of the bandage and pull to break the release liner. Remove release liner from bandage.

Option 1:



(3) Place the roll on the puncture site. Apply manual compression on the roll while the introducer is removed. To maintain manual compression on the roll, Step 4 may require assistance from a second health care provider.


(4) **DO NOT MOVE OR LIFT THE ROLL.** Center the adhesive bandage over the roll. Maintain manual compression on the roll. Stretch bandage around wrist until the adhesive edges overlap, adhering the tan colored tab to the patient's skin and overlapping the white tab to adhere on top of the tan tab. Note: Arterial punctures using large dilators will require longer manual compression time. Proceed to Step 5.

Option 2:


(3) Center the adhesive bandage over the puncture site. Stretch bandage around wrist until the adhesive edges overlap, adhering the tan colored tab to the patient's skin and overlapping the white tab to adhere on top of the tan tab.


(4) Insert the roll under the bandage at the puncture site, applying manual compression on the roll while the introducer is removed. Note: Arterial punctures using large dilators will require longer manual compression time. Proceed to Step 5.

(5) The adhesive bandage should be left in place for at least 30 minutes. Compression on the roll should then be gradually released by cutting every other elastic band on both sides over the next 30 minutes.

(6) The adhesive bandage's pressure should be completely released at one hour after application. To complete the pressure release, cut all the remaining elastic bands.

(7) The hemostatic roll may be left in place for up to 24 hours. NOTE: Health care providers are encouraged to continue to use the standard of care at their Institution. The standard of care includes manual compression at the puncture site, time to ambulation, and time to patient discharge.

PRECLINICAL DATA

Animal testing has been performed to demonstrate compatibility of the QuikClot and drug/induced anticoagulant therapy.

Study Name/Description	Pre-clinical evaluation of QuikClot in an anticoagulated swine model.																																
Objective	To evaluate the efficacy of kaolin, the active ingredient in Z-Medica's QuikClot in controlling bleeding in animals treated with common drugs such as Coumadin and Plavix.																																
Number of sites/investigators (OUS/US)	One site (US) and one investigator.																																
Number of animals/wounds	Ten pigs were divided in 2 groups (five animals each): One group received Plavix at 75 mg/day orally for > 5 days. The second group received daily doses of Coumadin until INR was > 3. 187 intra-abdominal vascular injuries (splenic, liver and mesenteric) were tested. Injuries consisted of surgically inflicted wounds with mixed arterial and venous bleeding. The size of the wounds were 5 cm in length and 3-5 mm in depth for liver and spleen and 2 mm of depth for the mesentery.																																
Procedure	In combination with manual pressure, QuikClot was applied to a series of wounds to the liver, spleen, mesentery and femoral artery and compared to standard surgical gauze. The study was prospective, open label, randomized 1:1 between QuikClot and control. Manual compression was applied over the wounds for 5 minutes and then released.																																
Success/Failure criteria	Success = complete bleeding cessation within 5 minutes. Failure = persistent bleeding at 5 minutes.																																
Study results	<p>Animals treated with Coumadin and QuikClot showed successful control of bleeding in 94.5% of wounds. Animals treated with Coumadin and control gauze showed control of bleeding in 24% of wounds. Animals treated with Plavix and QuikClot showed successful control of bleeding in 91.2% of wounds. Animals treated with Plavix and control gauze showed control of bleeding in 29.7% of wounds.</p> <p>QuikClot versus Control in Coumadin treated Pigs (n=5)</p> <table border="1"> <thead> <tr> <th></th> <th>Pass</th> <th>Fail</th> <th>Total</th> </tr> </thead> <tbody> <tr> <td>Test</td> <td>52</td> <td>3</td> <td>55</td> </tr> <tr> <td>Control</td> <td>9</td> <td>29</td> <td>38</td> </tr> <tr> <td></td> <td></td> <td></td> <td>93</td> </tr> </tbody> </table> <p>QuikClot versus Control in Plavix treated Pigs (n=5)</p> <table border="1"> <thead> <tr> <th></th> <th>Pass</th> <th>Fail</th> <th>Total</th> </tr> </thead> <tbody> <tr> <td>Test</td> <td>52</td> <td>5</td> <td>57</td> </tr> <tr> <td>Control</td> <td>11</td> <td>26</td> <td>37</td> </tr> <tr> <td></td> <td></td> <td></td> <td>94</td> </tr> </tbody> </table>		Pass	Fail	Total	Test	52	3	55	Control	9	29	38				93		Pass	Fail	Total	Test	52	5	57	Control	11	26	37				94
	Pass	Fail	Total																														
Test	52	3	55																														
Control	9	29	38																														
			93																														
	Pass	Fail	Total																														
Test	52	5	57																														
Control	11	26	37																														
			94																														
Adverse events	No animals died because of bleeding. The wounds where QuikClot did not control bleeding completely within 5 minutes required some additional manual compression time to achieve hemostasis.																																

CLINICAL DATA

Several peer reviewed clinical publications support the use of QuikClot in human patients on drug/induced anticoagulant therapy.

Study Name/Description	Trabattoni D, et al. A New Kaolin-based Hemostatic Bandage Use after Coronary Diagnostic and Interventional Procedures. Int J Cardiol. 2012;156(1):53-54.														
Number of sites/investigators (OUS/US)	One site (OUS) and three investigators.														
Number of subjects	40 subjects														
Inclusion/Exclusion	Patients undergoing diagnostic angiography or Percutaneous Coronary Intervention (PCI) via a femoral artery approach. Exclusion not described in publication.														
Procedure	Prospective single arm pilot trial of QuikClot in cardiac catheterization. Introducer sheath size 6 Fr (90%) or 7 Fr (10%). Femoral artery sheath removed once the ACT < 180 seconds. All patients were treated with QuikClot after arterial sheath removal.														
Study endpoints and assessment protocol	Complete bleeding cessation at 5 minutes and safe ambulation at 4 hours.														
Duration of follow-up	30 days														
Patient demographics	75% male, age 68+/-11 years														
Patient condition	<p>Patient undergoing diagnostic angiogram (62%) or PCI (38%) via femoral artery approach using 6F (90%) or 7F (10%) sheath.</p> <table border="1"> <thead> <tr> <th></th> <th>QuikClot (n=40)</th> </tr> </thead> <tbody> <tr> <td>LMW Heparin</td> <td>2.5%</td> </tr> <tr> <td>Aspirin + Clopidogrel</td> <td>27.5%</td> </tr> <tr> <td>Aspirin</td> <td>60%</td> </tr> <tr> <td>Aspirin + Warfarin</td> <td>5%</td> </tr> <tr> <td>IV Heparin</td> <td>38%</td> </tr> <tr> <td>No anticoagulation</td> <td>5%</td> </tr> </tbody> </table>		QuikClot (n=40)	LMW Heparin	2.5%	Aspirin + Clopidogrel	27.5%	Aspirin	60%	Aspirin + Warfarin	5%	IV Heparin	38%	No anticoagulation	5%
	QuikClot (n=40)														
LMW Heparin	2.5%														
Aspirin + Clopidogrel	27.5%														
Aspirin	60%														
Aspirin + Warfarin	5%														
IV Heparin	38%														
No anticoagulation	5%														
Study results	<p>Mean ACT value at hemostasis 138 ± 24 sec (range 95-186 sec). Mean cumulative hemostasis time 4.9±1.05 min. Diagnostic procedures 4.2±0.9 min. Interventional procedures 5.3±0.95 min. Ambulation time 4 h for all patients.</p>														
Adverse events	One PCI patient required extra compression time to achieve hemostasis and developed a small (< 5cm) hematoma.														

CLINICAL DATA

Study Name/Description	Politi L, et al. Randomized Clinical Trial on Short-Time Compression with Kaolin Filled Pad: A New Strategy to Avoid Early Bleeding and Subacute Radial Artery Occlusion after Percutaneous Coronary Intervention. J Internat Card. 2011; 24(1):65-72.																														
Number of sites/investigators (OUS/US)	One site (OUS) and ten investigators.																														
Number of subjects	120 subjects divided in 3 groups: Group 1 (QuikClot 15 minute compression, n=50), Group 2 (control – 15 minute compression, n=20), Group 3 (control 2 hour compression, n=50).																														
Inclusion/Exclusion	All patients undergoing transradial elective diagnostic or interventional coronary procedures between November 1, 2009 and January 31, 2010. Abnormal Allen's test before puncture excluded. Failure to provide written informed consent excluded.																														
Procedure	QuikClot was applied to the radial artery over the sheath which was then removed. Pressure was maintained for 15 minutes and then completely relaxed.																														
Study endpoints and assessment protocol	The main end-point was subacute Radial Artery Occlusion (RAO). The secondary end-point was failure of the closure technique (death, MI or major bleeding occurring in hospital). Groups 1 and 2: 15 minute assessment for bleeding. Group 3: 2 hour assessment for bleeding. All groups: Radial artery patency assessed at 24 hours using Barbeau's Test.																														
Duration of follow-up	Until patient discharge and follow-up visit. Follow-up at 6 months was done for patients who developed RAO.																														
Patient demographics	Age (years) Group 1 = 64.16 ± 11.53, Group 2 = 61.30 ± 14.22, Group 3 = 59.72 ± 14.23 Male Group 1 = 37 (74%), Group 2 = 14 (70%), Group 3 = 36 (72%) Weight (kg) Group 1 = 76.42 ± 11.13, Group 2 = 80.00 ± 14.96, Group 3 = 82.02 ± 13.26																														
Patient condition	<table border="1"> <thead> <tr> <th></th> <th>QuikClot</th> <th>Control Group 1</th> <th>Control Group 2</th> </tr> </thead> <tbody> <tr> <td>Aspirin</td> <td>46 (92%)</td> <td>19 (95%)</td> <td>50 (100%)</td> </tr> <tr> <td>Clopidogrel</td> <td>10 (20%)</td> <td>2 (10%)</td> <td>11 (22%)</td> </tr> <tr> <td>LMW Heparin</td> <td>8 (16%)</td> <td>5 (25%)</td> <td>5 (10%)</td> </tr> <tr> <td>Warfarin</td> <td>6 (12%)</td> <td>5 (25%)</td> <td>5 (10%)</td> </tr> <tr> <td>IV heparin</td> <td>100%</td> <td>100%</td> <td>100%</td> </tr> <tr> <td>No anticoagulation</td> <td>0%</td> <td>0%</td> <td>0%</td> </tr> </tbody> </table> <p>Several patients received multiple therapies</p>				QuikClot	Control Group 1	Control Group 2	Aspirin	46 (92%)	19 (95%)	50 (100%)	Clopidogrel	10 (20%)	2 (10%)	11 (22%)	LMW Heparin	8 (16%)	5 (25%)	5 (10%)	Warfarin	6 (12%)	5 (25%)	5 (10%)	IV heparin	100%	100%	100%	No anticoagulation	0%	0%	0%
	QuikClot	Control Group 1	Control Group 2																												
Aspirin	46 (92%)	19 (95%)	50 (100%)																												
Clopidogrel	10 (20%)	2 (10%)	11 (22%)																												
LMW Heparin	8 (16%)	5 (25%)	5 (10%)																												
Warfarin	6 (12%)	5 (25%)	5 (10%)																												
IV heparin	100%	100%	100%																												
No anticoagulation	0%	0%	0%																												
Study results	While a total of 150 patients were planned to be enrolled, the study stopped enrolling group 2 subjects after the 20th patient due to unethically high rates of bleeding. None of the patients enrolled in Group 1 (QuikClot) developed RAO the main outcome variable. Among patients enrolled in Group 2, RAO occurred in 1 case (5%) and among Group 3 in 5 cases (10%).																														
Adverse events	Active bleeding after compression removal: Group 1: 10 patients (20%) Group 2: 18 (90%) Group 3: 1 (2%) In all cases, hemostasis was achieved with a supplementary compression for 2 hours that did not produced any RAO in Group 1.																														

CLINICAL DATA

Study Name/Description	Trabattoni D, et al. A New Kaolin-based Haemostatic Bandage Compared with Manual Compression for Bleeding Control after Percutaneous Coronary Procedures; Eur Radiol, 2011 Aug 21(8): 1687-91.																						
Number of sites/investigators (OUS/US)	One site (OUS) and six investigators.																						
Number of subjects	200 subjects. Prospective randomized trial of QuikClot (n=100) vs. manual compression (n=100)																						
Inclusion/Exclusion	Undergoing angiography or Percutaneous Coronary Intervention (PCI) via a femoral approach. Patients with baseline INR > 1.4 excluded. Patients who had previous arterial access at the same femoral site within 30 days excluded.																						
Procedure	Femoral arterial sheath removed once the ACT < 180 seconds. Patients randomized to receive QuikClot or manual compression after femoral sheath removal. Patient ambulation at 4 hours.																						
Study endpoints and assessment protocol	Complete bleeding cessation at 5 minutes and safe ambulation at 4 hours.																						
Duration of follow-up	30 days																						
Patient demographics	Male 70% QuikClot vs. 60% control. Mean age (years) 65.7±13 QuikClot vs. 73.6±6.2 control. Weight (kg) 73.9±12 QuikClot vs. 71.2±15 control. Diagnostic procedure (n=98) or PCI (n=102). Introducer sheath size 6 Fr (90%) or 7 Fr (10%).																						
Patient condition	<table border="1"> <thead> <tr> <th></th> <th>QuikClot</th> <th>Control</th> </tr> </thead> <tbody> <tr> <td>LMW Heparin</td> <td>4%</td> <td>2%</td> </tr> <tr> <td>Aspirin + Clopidogrel</td> <td>29%</td> <td>26%</td> </tr> <tr> <td>Aspirin</td> <td>60%</td> <td>60%</td> </tr> <tr> <td>Aspirin + Warfarin</td> <td>7%</td> <td>3%</td> </tr> <tr> <td>IV Heparin</td> <td>51%</td> <td>49%</td> </tr> <tr> <td>No anticoagulation</td> <td>0%</td> <td>9%</td> </tr> </tbody> </table>			QuikClot	Control	LMW Heparin	4%	2%	Aspirin + Clopidogrel	29%	26%	Aspirin	60%	60%	Aspirin + Warfarin	7%	3%	IV Heparin	51%	49%	No anticoagulation	0%	9%
	QuikClot	Control																					
LMW Heparin	4%	2%																					
Aspirin + Clopidogrel	29%	26%																					
Aspirin	60%	60%																					
Aspirin + Warfarin	7%	3%																					
IV Heparin	51%	49%																					
No anticoagulation	0%	9%																					
Study results	Mean ACT value at hemostasis 146 ± 24 sec (range 98 – 198 sec). Hemostasis with QCI bandage 5.4 ± 1.5 min. Hemostasis with manual compression 25 ±15 min, p<0.001. No hemostasis failure in either group.																						
Adverse events	Major Bleeding: 1 patient in each group. Haematoma > 5 cm: 1 patient* QuikClot vs. 2 patients control. Pseudoaneurysm: 1 patient* QuikClot vs. 1 patient control. *Same patient																						

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